



When all you want is clear skin.



PicoWay™
The Clear Solution
For Pigmentation

Schedule a personal consultation
with your practitioner today!

- Fewer treatments
- Better clearance
- Less discomfort

SYNERON CANDELA®
www.syneron-candela.com

©2014. All rights reserved. PicoWay, Syneron, the Syneron logo are trademarks of Syneron Medical Ltd. and may be registered in certain jurisdictions. Candela is a registered trademark of the Candela Corporation. PB86721EN

SYNERON CANDELA®



I tried whitening creams, dermabrasion and other traditional laser treatments, but PicoWay was the only one to provide me with complete clearance with minimal side effects and discomfort.

DB, San Francisco



On a clear day....

Picoway® is the latest and most advanced laser solution for benign pigmentation such as freckles, age spots, lentigines, Nevus of Ota, Café-au-Lait and more. This breakthrough technology uses the shortest laser pulse on the market measured in picoseconds (trillionths of a second) to break down the pigment into the smallest particles.

Why PicoWay treatments?

- Effective on all pigmentation
- Clearance in fewer treatments
- Better results
- Less discomfort

How does the PicoWay work?

Laser light is targeted at unwanted pigment in ultrashort picosecond pulses. These pulses shatter the pigment into very tiny particles which are then easily eliminated by the body's natural processes. The PicoWay laser emits less heat energy than traditional lasers providing a safer, more comfortable treatment with fewer side effects.

Does the PicoWay treatment hurt?

The degree of discomfort varies with each person depending on the nature and location of the pigmented lesions on the body.

What should you expect and what precautions should be taken throughout this process?

During your first visit, your skin's reaction to the laser will be tested to determine the most effective treatment.

Immediately following treatment, cold or cool compresses may be applied to soothe the treated area. A topical antibiotic ointment may be recommended. In the first few days, take care to avoid scrubbing the area or using abrasive skin cleansers.

How many PicoWay treatments will I require?


The effectiveness of the treatment depends on many factors such as the location, depth, size and type of pigmentation. Some lesions will require several treatments for optimal results, where some can be resolved in a single treatment.

Are PicoWay treatments safe, and are there any side effects?


PicoWay treatments selectively treat the melanin pigmentation while leaving surrounding skin undamaged. After treatments, patients will experience temporary discoloration or reddening of the skin around the treatment site. Some patients may also experience a temporary lightening or darkening of the skin for a few months following treatment. The skin will typically resolve to normal pigmentation color.

Can the PicoWay also treat tattoos and other discolorations?

Yes, the PicoWay treatment is also safe and effective for treatment of multicolored tattoos.



a. Pigmentation before PicoWay treatment



b. PicoWay technology shatters the pigment into minuscule particles that can easily be removed by the body.