

TATTOO REMOVAL

PICOWAY LASER TATTOO REMOVAL

PICOWAY LASER BY CANDELA MEDICAL

The most advanced picosecond laser used to treat tattoos of all colours. The PicoWay shatters tattoo ink into miniscule particles which is cleared by the body's immune system. It's safe for darker skin types and with less discomfort during and after treatment, tattoo removal sessions can be spaced closer together and your tattoo is treated effectively and in fewer treatments.

Why is PicoWay superior to other laser tattoo removal?

PicoWay has the shortest pulses and the highest peak power of any tattoo removal laser. This means better ink clearance (all colours) and less heat is delivered to the skin so there's less discomfort and faster healing.

For best results, treatments are spaced 6-8 weeks apart.

FREQUENTLY ASKED QUESTIONS

How many treatments will I need?

The effectiveness of laser tattoo removal is dependent on several variables such as age of the tattoo, location, depth, size, the type of ink used among other factors. Consultations with our experienced therapists are complimentary. They will assess your tattoo and advise the number of treatments required for fading or complete removal.

What if I only want to fade my tattoo?

Opt for one of our 'cover up' packages if you would like to refresh the ink in your tattoo or 'cover up' an existing tattoo.

What if I need more treatments?

Don't pay full price. Prepay your recommended package and if you need additional treatments, you can pay as you go and save 10% on PAYG prices.

SIZE GUIDE



XS

Up to size
of a postage
stamp



S

Up to size of
a business
card



M

Up to size of
\$5 note



L

A6,
¼ sleeve



XL

A5,
½ sleeve

PRICING

PACKAGE PRICING

Save up to 30%, no more to pay!

TREATMENTS

PAYG
per session

4
cover up

6
popular

8

10

CONSULTATION

complimentary

EXTRA SMALL (XS)

130

440

625

780

910

SMALL (S)

230

780

1100

1380

1600

MEDIUM (M)

330

1120

1580

1980

2310

LARGE (L)

430

1460

2060

2580

3010

EXTRA LARGE (XL)

630

2140

3020

3780

4400